

# Intech

## Bottle Top Dispenser



Intech is proud to offer its NEW BOTTLE TOP DISPENSERS which is an unique combination of affordable pricing and high performance.

Simple design clubbed with variety of useful features makes Intech Bottle top Dispensers an obvious choice for all routine dispensing needs.

High precision and accuracy is ensured through several stages of strict quality checks during manufacturing process. Each instrument is individually calibrated in accordance with ISO 8655 standards and comes with Individual Calibration certificate.

### SALIENT FEATURES

- Three models with adjustable volumes :  
(a) 0.25 - 2.5 ml (b) 0.5 - 5 ml (c) 1 - 10 ml
- Easy to use, simple construction and a good chemical compatibility.
- Each instrument comes with three adapters to fit in most of the laboratory reagent bottles. Standard threads : 30mm, Adapters : 25mm, 36mm and 45mm
- Smooth, effortless plunger movement.
- Safety discharge system reduces the risk of accidental dispensing.
- Nozzle cap prevents any unwanted drop on the work space.
- Instrument is fully autoclavable at 121°C , 15psi. No Need for dis-assembling the instrument for autoclaving.
- Easy to dis-assemble for cleaning and servicing.
- In-lab, user re-calibration is possible in compliance with GLO/ISO Norms.

### PRODUCT SPECIFICATIONS

Piston	Polypropylene
Cylinder	Glass
Housing	Polypropylene
Volume adjustment	Locking rocker switch
Valve assembly	Glass Ball & Seat
Discharge Assembly	Teflon
Delivery tube	Polypropylene
Calibration	Individually calibrated and certified
Accuracy & Reproducibility	In accordance with ISO 8655 standards
Compatibility	Excellent- Suitable for most reagents except HF
Bottle Thread Size	30mm
Additional Adaptors	25mm, 36mm, 45mm

\* Bottle top dispenser available in four colours.

\* Design & Specifications are subject to change without any prior notice.

### SPECIFICATIONS & ORDERING INFORMATION

#### Bottle Top Dispenser

Model No.	Vol.		Accuracy		CV	
	Range	Increment	±%	±%ml	±%	±%ml
IN-D1	0.25 - 2.5 ml	0.05 ml	0.6	0.012	0.2	0.004
IN-D2	0.5 - 5 ml	0.1 ml	0.6	0.030	0.2	0.010
IN-D3	1 - 10 ml	0.2 ml	0.6	0.060	0.2	0.020

Perfection in  
Laboratory Science

Telefax : 0091-172-2563873, 0091-9815611479  
international@rediffmail.com, exports@micropipettes.net  
www.micropipettes.net

An ISO 9001:2000 Certified Company

## Bottle Top Dispenser



# Intech

An ISO 9001:2000 Certified Company